

第四届苏州国际医学影像研讨会

Deep Vision for Human Vision

Yu Qiao¹ and Zhe Wang²

1. Shenzhen Institutes of Advanced Technologies, Chinese Academy of Sciences

2. Sensetime Limited









SENSETIMIC 商汤科技

SenseTime: Who We Are

SenseTime is an AI company focusing on:

① Computer Vision

2 Deep Learning



SenseTime developed a proprietary deep-learning platform, the only one in China

1

SENSETIME 商 汤 科 技

SenseTime – The Largest Al Unicorn in China



Leader in R&D, Talent, and Business

Research Achievements

- Core team has 20+ years of research experience. One joint lab was named 2016 Top 10 AI Lab in the World
- 400+ research papers published at top conferences
- Beat Facebook to be the first team that achieved better accuracy than human eye with facial recognition AI
- Competed with Google at ImageNet Large Scale Visual Recognition Challenge

CVPR2016

kaggle

IM AGENET

3

The Only Deep Learning Platform in China

- Developed deep learning platform including both software framework and supercomputing system
- Significantly reduced the R&D cost and shortened the time needed in training highly sophisticated models



Deep Talent Pool

- The largest deep learning research team in Asia
- 18 professors, 120+ PhDs, and 300+ masters/bachelors from top universities in the world
- One of the largest AI research teams in the world



#1 in AI Commercialization

SenseTime's 2016 revenue topped the industry. Customers:

- Video surveillance providers
- Banks and financial institutions
 - Smart phone OEMs
- Mobile application developers
 - Internet service providers
 - Robot makers
- Government agencies & public security authorities

SENSETINIC 商 汤 科 技

Computer Vision v.s. Human Vision

V. S.



Human vision

- Evolution result
- (~520m yrs)
- Bio-chemical events
- Powerful, robust and ... reliable



Computer vision

- A sub area of AI
- (~60 yrs)
- Electrical computation
- Still needs requirements

Computer vision vs human vision





Three basic problems of computer vision

SENSETINC 商 汤 科 技

Strengths of deep learning

Big data

Powerful computation

Deep models



Surpass human performance

7





ImageNet: Object Classification DL 97.21%: Human 94.9%

LFW: Face Verification DL 99.6%: Human 97.3%

Can deep learning help, enhance or assistant human vision?



SENSETIME 商 汤 科 技

Deep Learning for Human Vision



Super-resolution & Image enhancement



Diabetic Retinopathy Diagnosis



Glaucoma Diagnosis



Skin lesion analysis

Super-resolution and Image Enhancement



Image enhancement for low-light photos



SENSETINIC 商汤科技

Image Super-resolution







Image deblurring





Deblurring for medical image





SENSETINIC 商 汤 科 技

Deblurring for medical image





Paper and Challenge Results

SENSETIME 商汤科技

SENSETIME

- "Deep Convolutional Neural Network for Image Deconvolution", NIPS 2014
- "On Vectorization of Deep Convolutional Neural Networks for Vision Tasks", AAAI 2015
 - "Deep Edge-Aware Filters" , ICML 2015
 - "Shepard Convolutional Neural Networks", NIPS 2015



- "Learning a Deep Convolutional Network for Image Super-Resolution", ECCV 2014 (539 citations)
- "Image Super-Resolution Using Deep Convolutional Networks", TPAMI 2015 (461 citations)

ICDAR





(qB)

IAPR (2)

2nd in CVPR 2017 NTIRE Image supersolution challenge

1st in ICDAR 2015 Blurring Text Classification challenge

SENSETIME 商汤科技

Diabetic retinopathy Diagnosis

Diabetic Retinopathy is the leading cause of blindness in the working age population in developed countries. In 2015, it is estimated to affect over 93million people.





3:Severe

4:Proliferative

```
2:Moderate
```

SENSETIME

科

技

汤

Google has published a paper on JAMA in 2017. They asked 54 ophthalmologists to label a total of 130 thousand images, and train a neural network with this dataset, which achieved a sensitivity of 90.3% and specificity of 98%.

We proposed Zoom-in-Net achieved state-of-the-art performance on DR classification. Trained with only image-level labels, it can also accurate localize lesion regions with a recall score of 80% with 4 proposals per image.

Observation

• The size of a typical retinal image of fundus is about 4000*4000 ,while the size of a typical input image of a modern CNN is smaller than 500*500.

SENSETIME

朻

顸

汤

苗

• The input size is limited by the GPU memory as well as the network complexity.



To get a sense of what the size difference is

Our Solution is Zoom-in-Net, which mimics the process of clinician' s behavior:

SENSETIME

科

汤

商

技

- 1. Scan the whole image to identify the suspicious regions
- 2. Automatically zoom in these regions to verify



Zoom in suspicious regions for details:

Sampling highest response regions on the gated attention maps in a greedy manner.



From left to right: image, gated attention maps of level 1-4 and the selected regions of the image.

Performance on Kaggle challenge and Messidor dataset

Algorithms	val set	test set	Method	AUC	Acc.	Method	AUC	Acc.
Min-pooling [45]	0.86	0.849	Lesion-based [50]	0.760	-	Splat feature/kNN [52]	0.870	-
o_O [46]	0.854	0.844	Fisher Vector [50]	0.863	-	VNXK/LGI [48]	0.870	0.871
Reformed Gamblers [47]	0.851	0.839	VNXK/LGI [48]	0.887	0.893	CKML Net/LGI [48]	0.862	0.858
M-Net	0.832	0.825	CKML Net/LGI [48]	0.891	0.897	Comprehensive CAD [51]	0.876	-
M-Net+A-Net	0.837	0.832	Export A [51]	0.91	-	Expert A [51]	0.922	-
Zoom-in-Net	0.857	0.849	Expert R [51] Expert B [51]	0.94 0.92	-	Expert B [51]	0.865	-
Ensembles	0.865	0.854	Zoom-in-Net	0.957	0.911	Zoom-in-Net	0.921	0.905

 Table 2.1: Comparison to top-3 entries on Kaggle's challenge.

 Table 2.2: AUC for referable/nonreferable

 Table 2.3: AUC for normal/abnormal

Sansatillic

顸

汤



SENSETIME

汤科

商

技

- Zoom-in-Net v.s. general physician v.s. optometrists¹
 0.865 v.s. 0.838 v.s. 0.719
- 83% Recall for Zoom-in-Net with four proposals per image (left)
- Automatically discover meaningful lesion types (right)

SENSETINIC 商汤科技

Glaucoma Diagnosis

- Glaucoma diagnosis
- Sensetime collaborated with Shenzhen Institute of Advanced Technology and ten ophthalmic hospitals.

Sansatillic

技

汤科

商

- Over 7k visual fields, 10k OCT images and 10k fundus images, the largest dataset for comprehensive diagnosis of glaucoma.
- Based on the visual reports, we designed a deep-learning-based algorithm which outperforms human experts.



- Glaucoma diagnosis
- Our system can automatically verify the reliability of a visual field report by its fixation losses, False positive and False negative errors. And then it extracts the Pattern Deviation (PD) plot and process it with a convolution neural network.

Sansatime

汤科

商

技



• Network structure

• The backbone of the network is Based on VGG16. We modify the output dimension of the penultimate layer fc7 from 4096 to 200 so that it is more comparable to the dimension of the three parameters VFI/MD/PSD.

SENSETIME

汤 科

商

技

• To fuse the potential information from VFI/MD/PSD, the 200-dimension output and the three parameters are concatenated to form a fused vector, which is used for final prediction.



- Glaucoma diagnosis
- Left figure shows baseline statistics
- Right figure shows the performance on a test dataset of 300 PD plots, compared to traditional approaches and nine ophthalmologists.

	Non-glaucoma Group	Glaucoma Group	*P Value	
No. of images	1623	2389	_	
age	47.2 (17.4)	49.2 (16.3)	0.0022	
left/ri ght	635/919	607/911	0.6211	
VFI	0.917 (0.126)	0.847 (0.162)	0.0001	
MD	-5.0 (23.5)	-9.0 (44.8)	0.0039	
PSD	3.6 (3.3)	6.7 (22.2)	0.0001	

Mean and standard deviation are provided. *Comparison between normal group and glaucoma group (unpaired t test for numerical data and chi-square test for categorical data)



SENSETIME

科

汤

商

技

• Now we are developing **a cloud-based platform** for clinical trials. We hope this A.I. system can really help improve the diagnosis accuracy

• We are also working on integrating OCT images and fundus images in the system.

SENSETIME

技

汤

商

• Different modalities, better performance



Visual Field + Fundus Image + OCT

SENSETIME 商汤科技

Skin Lesion Analysis

Skin Lesion Analysis

- Skin lesion diagnosis
- With big data and deep learning, A.I. can perform as good as human experts¹.



SENSETIME

朻

汤

商

技

• Currently, our system outperforms the state-of-the-art methods on ISIC Skin challenge datasets on both the classification task and segmentation task.

1.Dermatologist-level classification of skin cancer with deep neural networks, Nature 2017

Skin Lesion Analysis

• Skin lesion classification and segmentation are two highly correlated tasks. However, their relationship is not fully utilized in previous methods.

SENSETIME

汤

- We designed a multi-task deep convolution neural network architecture to solve both tasks simultaneously.
- A feature passing module is proposed for message transmission between the two networks.



Fig. 1. Workflow of the proposed algorithm

• Currently, our system achieves the accuracy of 81% and a mIoU of 0.79 on the ISIC challenge (both best).

1. Dermatologist-level classification of skin cancer with deep neural networks, Nature $2017\,$

Skin Lesion Analysis

Method	AC(cls)	mAP(cls)	JA(seg)	AC(seg)	DI(seg)
Base	0.772	0.699	0.779	0.940	0.862
Multi-task	0.779	0.712	0.780	0.940	0.863
Ours	0.801	0.747	0.787	0.944	0.868

 Table 1. Performances compared with our base model and
 multi-task network without feature passing module.

Table 2. Classification performance compared with other Table 3. Segmentation performance compared with other state-of-art classification models.

Method	AC(cls)	mAP(cls)
AlexNet	0.775	0.706
VGG16	0.789	0.727
ResNet101	0.772	0.699
Ours	0.801	0.747

state-of-art segmentation models.

SENSETIME

商

汤科技

Method	JA(seg)	AC(seg)	DI(seg)
MtSinai	0.765	0.934	0.849
U-net	0.741	0.926	0.822
Deeplab-ResNet101	0.779	0.940	0.862
Ours	0.787	0.944	0.868

Conclusions

Deep learning for human vision in many ways
➢ Better view of photos (enhancement, super-resolution, …)
➢ Diagnosis of eye related disease (Glaucoma, Diabetic
Retinopathy,…)
➢ Scan huge image datasets (medical, surveillance, …)

SENSETIME 商汤科技 SENSETIME 商汤科技

THANK YOU!

Follow Sensetime on WeChat

